

שער יוקנעם בנין תבור
ת.ד. 595 יוקנעם עילית 20692.
טל. 04-9099100
פקס. 04-9591982
דואר אלקטרוני: Itsik.benatar@MEDiGAL.com

P.O.X 595
YOKNEAM ILLIT 20692
Tel. 972-4-9099100
Fax. 972-4-9591982
E-mail: Itsik.benatar@MEDiGAL.com

Preventive Maintenance Form- HC10

Customer Name: _____

Date: ____ / ____ / ____

Instrument s/n: _____

Next Maintenance: ____ / ____

Last Update: 20.3.2018

Comment	Value:	Checked/ Done:	Remarks
1. External checks			
1.1 General cleaning			
1.2 Monitor cleaning			
1.3 Tubes check (back)			
1.4 Tube Adaptor holder cleaning			ניקוי האדפטורים והשקע בדלתית
1.5 Run Decontamination Procedure			
2. Software			
2.1 Software version			Version: _____
2.2 Preform Self-Test			Maintenance > Diagnostics > Self-test > Start
2.4 Statistics (Service User), check if there are frequent problems			Maintenance > Diagnostics > Service > Diagnostics> Statistics > Statistics
2.5 from statistics, Number Of Successful Measurements			Maintenance > Diagnostics > Service > Diagnostics> Statistics > Statistics
2.6 Reset- איפוס רק לאחר אישור ממנהל השירות	Reset?		Maintenance > Diagnostics > Service > Diagnostics> Statistics > Statistics
3. Hardware			
3.1 Dilutor Piston- Checking & Lubricating Tips			
3.2 Washing Station- Cleaning & check			
3.3 Needle- Cleaning & check			
3.4 Test All Motors			Maintenance > Diagnostics > Service > Trouble-shooting> Test All
3.5 Run Decontamination Procedure (according to Samsung instructions)			
4. Functional checking's			
4.2 Run QC Batch (Do Calibration if necessary)	לצרף לדוח		

Service Engineer Name: _____

Service Engineer signature: _____

Decontamination Procedure

1. Execute a complete preparing for shipment procedure using distilled water when prompted. System will turn off.
2. Restart system.
3. Connect Hypo-cleaner (D-701) solution or a 10% dilution (1:9) of Clorox (bleach) to all reagent input lines except waste, and perform priming on all reagent inputs.
4. Run a Hard Cleaning procedure on the instrument. Follow the onscreen instructions.
5. Run 3 x blank measurements
6. Leave the instrument 30 minutes untouched
7. Execute a complete preparing for shipment procedure again using distilled water when prompted. System will turn off.
8. Restart system.
9. Replaced the complete diluent reagent line on the analyzer including the cap and pickup tube.
10. Connect all system reagents, making sure to use a new unopened container of diluent.
11. Prime all reagents on the system.
12. Execute blank measurements. (It is normal if the first few runs are elevated)

Service Engineer Name: : _____ **Service Engineer signature:** _____