	<h1>Work Procedure</h1>	WP-0278-15
		Page 1 / 4
Created by: SAV	Approved: AQ	11.11.2013

Procedure name

WP-0278 Sepax final test

Select the type of device

	article N°	Denomination	Serial Number
<input type="checkbox"/>	12000	Sepax S-100	
<input type="checkbox"/>	13000	Sepax RM	

Tools

Multimeter: _____ SecuTest: _____
 Test syringe 2ml with luer lock / Black filter / Test chamber / Empty test chamber / Kit with water/Torque wrench

Test Procedure

1. Customer parameters:

Check the **Service Menu** parameters

Serial Number: _____
 Driver printer: _____
 Piston position offset: _____

Toggle flags activated: Fan Monitoring: ☐
 Write summary data: ☐
 Enable Printer: ☐
 Covers Temp Sensor: ☐

Check the **Settings Menu** parameters:

AutoPrintData: ☐
 Date and time are correct: OK: ☐
 Language: _____

Check the **Traceability Setup** (if activated and in case of software update):

N/A ☐
 Enable Trace ID: ☐
 Manual ID input: ☐
 ISBT 128 only: ☐

Take the **Protocol** parameters in case of an upgrade of protocol or MAP:

Parameters: _____	Traceability ID: _____
_____	_____
_____	_____
_____	_____

2. Air tests

Attach the test syringe 2ml with luer lock to the line pressure detector.

+2ml of air _____	(1000 < value < 1344 mbar)	Line+ OK: <input type="checkbox"/>
-2ml of air _____	(440 < value < 660 mbar)	Line- OK: <input type="checkbox"/>

Install the empty test chamber in the centrifuge and connect to the line pressure detector.
 Run the air pump to create:

0 mbar _____ chamber _____ line	(ΔMax 20mbar)	ΔChr Line0 OK: <input type="checkbox"/>
-300 mbar _____ chamber _____ line	(ΔMax 20mbar)	ΔChr Line- OK: <input type="checkbox"/>
+500 mbar _____ chamber _____ line	(ΔMax 20mbar)	ΔChr Line+ OK: <input type="checkbox"/>

Apply a pressure of 500 mbar, close the stopcock on the line sensor and monitor the pressure drop for 1 min.

Starting pressure _____ mbar **Pressure after 1min** _____ line (ΔMax 20mbar) Line drop OK: ☐

Install the test chamber in the centrifuge.


Apply maximum pressure (ELv 4000) _____	(1500 < value < 1800 mbar)	Chr max OK: <input type="checkbox"/>
Apply maximum vacuum (ELv 4000) _____	(value < -700 mbar)	Chr min OK: <input type="checkbox"/>
Apply vacuum or pressure (ELv 0) _____	(50 < value < 150 mbar)	Chr 0 OK: <input type="checkbox"/>

Check if nothing touches the compressor and create any vibrations!!

Perform the test below only in case of sealing system exchange. Otherwise choose N/A

Apply a pressure of 500mbar and block the circuit. Monitor the pressure drop during 1 minute

Starting pressure _____ mbar **Pressure after 1 min** _____ mbar (ΔMax 35mbar) Chr drop OK: ☐ N/A : ☐

	<h1>Work Procedure</h1>	<h1>WP-0278-15</h1>
Created by: SAV		Page 2 / 4 11.11.2013

Install the test chamber in the centrifuge, apply a pressure of 500mbar, block the pressure and run the centrifuge at 4000rpm.

Wait 2 min and then monitor the pressure drop during 1 min.

Pressure at 2 min _____ mbar **Pressure at 3rd min** _____ Chr (Δ Max 30mbar) Chr drop OK: ☐

3. Chamber volume detector tests

Adjust the ref 2.5V for piston captor with 2 CDD

Empty centrifuge: _____ Pixels (=1) OK: ☐

Install the test chamber in the centrifuge.

PNC Position 1: _____ Pixels (69 < value < 76) OK: ☐

PNC Position 2: _____ Pixels Δ of position 1-2 (390<value<410) _____ Pos 1-2 OK: ☐

PNC Position 3: _____ Pixels Δ of position 2-3 (390<value<410) _____ Pos 2-3 OK: ☐

Put the value of position 1 in the Service menu at the position piston offset choice OK: ☐

Go back to the man_sepax to verify if the value is equal to 65 at the P value OK: ☐

4. Optical line sensor tests

Switch ON the Sepax and wait 10 min.

Without filter

Red LED _____ V (5.5 < value < 6.0V) (5.0 < value < 5,7V) Red WF OK ☐

Green LED _____ V (5.5 < value < 6.0V) (5.0 < value < 5,7V) Green WF OK ☐

Black filter

Red LED _____ V (0.0 < value < 0.2V) (0.0 < value < 0.2V) Red WF OK ☐

Green LED _____ V (0.0 < value < 0.2V) (0.0 < value < 0.2V) Green WF OK ☐

5. Security tests

Check the splash detector OK: ☐

Check the three chamber cover detectors OK: ☐

Check the two temperature sensors _____ (0.4 < C < 0.7) _____ (18°C < I < 40°C) OK: ☐

For the cover temperature sensor, pass the finger on the reading window. The value must change OK: ☐

6. Electrical power adjust

Check different voltage on the top of each fuse.

GND

+ 5.1V Fuse 1 _____ (5.0 < value < 5.2) DC1 OK: ☐

+ 15.1V Fuse 2 _____ (14.8 < value < 15.3) DC2 OK: ☐

-15V.1 With Lambda power supply and -12.1 With XP-power supply or Sepax with SN> 653:

- 15.1V _____ (-15.3 < value < -14.8)

- 12.1V Fuse 3 _____ (-12.3 < value < -11.8) DC3 OK: ☐

+ 26.1V Fuse 4 _____ (25.5 < value < 26.3) DC4 OK: ☐

+ 48.1V Fuse 5 _____ (47.5 < value < 48.3) DC5 OK: ☐

7. User function tests

There must be **three beeps** at start-up OK: ☐

Check that all the **LCD characters** are displayed OK: ☐

Note the **software** version n° : _____

Test the **touch panel** _____ OK: ☐

Switch on the **blue LEDs** OK: ☐

Main fan and power fan turn OK: ☐

Write the switch configuration on the interface board.


(SW2 OFF = ETEL motor without SepaxNet (COM1) / SW2 ON = Etel motor with SepaxNet (COM4))

(SW3 OFF = old piston captor 2CCD / SW3 ON = new piston captor 1CCD)

(SW4 OFF = ETEL motor / SW4 ON = SERVIDA motor)

(SW1 always OFF / SW5 always ON)

	OFF	ON
SW1	<input checked="" type="checkbox"/>	<input type="checkbox"/>
SW2	<input type="checkbox"/>	<input type="checkbox"/>
SW3	<input type="checkbox"/>	<input type="checkbox"/>
SW4	<input type="checkbox"/>	<input type="checkbox"/>
SW5	<input type="checkbox"/>	<input checked="" type="checkbox"/>

	Work Procedure	WP-0278-15
		Page 3 / 4 11.11.2013
Created by: SAV	Approved: AQ	

8. Centrifuge Tower Test – UCB protocol compatibility

Centrifuge tower alignment . The centrifuge should be perfectly centred	OK:	<input type="checkbox"/>
Centrifuge feet stability	OK:	<input type="checkbox"/>
Centrifuge upper deck (with plastic ring) is not touching (or not too close) to the device enclosure	OK:	<input type="checkbox"/>
The distance must be minimum 1mm to the U panel	OK:	<input type="checkbox"/>
Check if the centrifuge covers are neither broken nor bent, and do not create any vibrations	OK:	<input type="checkbox"/>
No internal component touching each other's (especially compressor, air buffer and centrifuge tower)	OK:	<input type="checkbox"/>

9. Stepper motor tests

Run the protocol "test stpck pos " with stop cock ramp and choose time 60s	OK:	<input type="checkbox"/>
With the torque wrench , control the stopcock blocked in the both direction	OK:	<input type="checkbox"/>

Electrical test according to IEC 62353

10. Electrical earth tests

Protective Earth Resistance / RCP mode

With the Secutest check the earth connection at 200mA. Test the following point. Max value <0.3Ω

The bolt under the rear power connection	OK:	<input type="checkbox"/>
The rear fan grill	OK:	<input type="checkbox"/>
A handle fixing bolt	OK:	<input type="checkbox"/>
The centrifuge earth bolt	OK:	<input type="checkbox"/>

11. Insulation test

Insulation Resistance / R-ISO mode

Connect the Sepax to the Secutest with 500V. Test the bolt of the centrifuge cover. Min Value >70M Ω

Centrifuge cover's bolt	OK:	<input type="checkbox"/>
-------------------------	-----	--------------------------

12. Earth leakage test

Leakage / ICD mode

Current Normal Condition / Differential	_____	Max value < 0.5mA	OK:	<input type="checkbox"/>
Single Fault Current / Alternative Equ.	_____	Max value < 1 mA	OK:	<input type="checkbox"/>

13. Current and voltage tests

In function mode

Current Consumption	_____	Max value 0.65A Max value 1A	@230Vac @110Vac	OK:	<input type="checkbox"/>
---------------------	-------	---------------------------------	--------------------	-----	--------------------------


Other tests

14. Protocol tests

Run the protocol test service using a kit with water	OK:	<input type="checkbox"/>
Attach the test file printout (if printer available)		

15. Visual control

Check the exterior appearance of the Sepax	OK:	<input type="checkbox"/>
---	-----	--------------------------

	<h2>Work Procedure</h2>	WP-0278-15
		Page 4 / 4
Created by: SAV	Approved: AQ	11.11.2013

16. AS-610 Traceability

Installed	YES: <input type="checkbox"/> NO: <input type="checkbox"/>
Check the functionality of the barcode reader reading a barcode test (see below)	OK: <input type="checkbox"/> N/A: <input type="checkbox"/>
Check the thermal printer exterior appearance and functionality	OK: <input type="checkbox"/> N/A: <input type="checkbox"/>
Check the printer power supply and power cable exterior appearance	OK: <input type="checkbox"/> N/A: <input type="checkbox"/>
Check the Sepax Net printing	OK: <input type="checkbox"/> N/A: <input type="checkbox"/>



17. Comments

18. Parts exchanged

Article Number	Designation	Old SN	New SN

In case of battery exchange	Reference	Manufacturing date

Tests conclusion

Tests passed for the Sepax

PASS: ☐
 NOT PASS: ☐

Approval

Date: _____ Signed: _____