	Work Procedure Maintenance	IT-1679-10
		Page 1 / 4 21.01.2014
Created by: SAV	Approved: AQ	

Procedure name

IT-1679 Sepax 2 Work Instructions

Select the type of device

	article N°	Denomination	Serial Number
<input type="checkbox"/>	14000	Sepax S-100 (2)	
<input type="checkbox"/>	14100	Sepax RM (2)	

Tools

SecuTest SN: _____

2ml syringe with luer lock / Pressure test chamber / Empty test chamber / Torque wrench / Water kit test

Test Procedure GMAP software

1. User function tests

Test the Stop Button (Hardware)	OK:	<input type="checkbox"/>
Check that the Touch screen displays well, touch-screen function works, backlight	OK:	<input type="checkbox"/>
Check the accuracy of the touchscreen display (calibrate it if needed)	OK:	<input type="checkbox"/>
Check that the blue LEDs are switched on	OK:	<input type="checkbox"/>

2. Customer parameters

Check the **Service Menu** Parameters

GMAP version	_____
GMAP build version	_____
Windows version	_____
Serial number	_____
Printer driver	_____
Piston position offset	_____
Debug => Send Logfile	<input type="checkbox"/>
Patfile Unicode	<input type="checkbox"/>
Page format	_____
SepaxNet host (If applicable)	_____

Toggle flag activation

Fan monitoring	<input type="checkbox"/>
Cover temp. sensor	<input type="checkbox"/>
Show menu button	<input type="checkbox"/>

Check the **Settings Menu** Parameters

Auto print data	<input type="checkbox"/>
Date and time are correct	<input type="checkbox"/>
Language:	_____
Sound volume:	_____
Filename donation ID	<input type="checkbox"/>


Check the **Traceability Setup** (if activated and in case of software update):

N/A	<input type="checkbox"/>
Enable trace ID	<input type="checkbox"/>
Force Input ID	<input type="checkbox"/>
ISBT128 only	<input type="checkbox"/>

In case of update/reinstallation of protocol, note the **protocol's parameters** below

Parameters: _____

Traceability ID: _____

	Work Procedure Maintenance	IT-1679-10
Created by: SAV		Page 2 / 4 21.01.2014

Test Procedure GMAN software or Automated Device Inspection Protocol

Automated Device Inspection Protocol V140

If Automated Device Inspection Protocol is used, skip chapters 3 to 8. Printout report must be then attached to this working procedure.

3. Electrical power adjust

Check the following output voltages in the **Manual software** Sepax 2

+ 48.00V	Fuse 1 (4A)	_____	(47.4 < value < 48.6)	DC1 OK:	<input type="checkbox"/>
+ 24.20V	Fuse 2 (2,5A)	_____	(23.7 < value < 24.7)	DC2 OK:	<input type="checkbox"/>
- 12.05V	Fuse 3 (1A)	_____	(-11.55 > value > -12.55)	DC3 OK:	<input type="checkbox"/>
+ 12.05V	Fuse 4 (250mA)	_____	(11.55 < value < 12.55)	DC4 OK:	<input type="checkbox"/>
+ 5.05V	Fuse 5 (4A)	_____	(4.8 < value < 5.3)	DC5 OK:	<input type="checkbox"/>

"If the values is showing a short peak as "minus" (digit over flow), then it can be ignored"

4. Air tests

Attach the test **2ml syringe with luer lock** to the line pressure detector.

+2ml of air	_____	(1200 < value < 1500 mbar)	Line+ OK:	<input type="checkbox"/>
-2ml of air	_____	(-650 < value < -450 mbar)	Line- OK:	<input type="checkbox"/>

Install the Empty test chamber in the centrifuge and connect the stopcock to the line pressure sensor.

Run the air pump to create: (When PRS/VAC required values are reached, use the BLK button to block the circuit).

0 mbar	_____	chamber	_____	line	(ΔMax 20 mbar) ΔChr Line OK:	<input type="checkbox"/>
-300 mbar	_____	chamber	_____	line	(ΔMax 20 mbar) ΔChr Line- OK:	<input type="checkbox"/>
+500 mbar	_____	chamber	_____	line	(ΔMax 20 mbar) ΔChr Line+ OK:	<input type="checkbox"/>

Apply a pressure of **~500 mbar**, close the stopcock and monitor the pressure drop for **1 min**.

Starting pressure	_____	mbar	Pressure after 1min	_____	line (ΔMax 5mbar)	Line drop OK:	<input type="checkbox"/>
-------------------	-------	------	---------------------	-------	-------------------	---------------	--------------------------

Install the Pressure test chamber in the centrifuge.

Apply maximum pressure	(ELv 4000)	_____	(1500 < value < 1900 mbar)	Chr max OK:	<input type="checkbox"/>
Apply maximum vacuum	(ELv 4000)	_____	(value < -650 mbar)	Chr min OK:	<input type="checkbox"/>
Apply vacuum or pressure	(ELv 0)	_____	(50 < value < 150 mbar)	Chr 0 OK:	<input type="checkbox"/>

Perform the test below only in case of sealing system exchange. Otherwise choose N/A

Apply a pressure of 500 mbar and block the circuit. Monitor the pressure drop during 1 minute

Starting pressure	_____	mbar	Pressure after 1 min	_____	mbar (ΔMax 35 mbar)	Chr drop	OK:	<input type="checkbox"/>	N/A:	<input type="checkbox"/>
-------------------	-------	------	----------------------	-------	---------------------	----------	-----	--------------------------	------	--------------------------

With the Pressure test chamber inside the centrifuge apply a pressure of **~500 mbar** and run the centrifuge at **4000 rpm**. Wait **2 min** and then monitor the pressure drop during **1 min**.

Pressure at 2 nd min	_____	mbar	Pressure at 3 rd min	_____	mbar (ΔMax 30 mbar)	Chr drop OK:	<input type="checkbox"/>
---------------------------------	-------	------	---------------------------------	-------	---------------------	--------------	--------------------------

5. Chamber volume sensor tests

Check the not Corrected value (**NC**) in the Manual Sepax without chamber

Offset Empty centrifuge :	_____	Pixels	(≤10)	OK:	<input type="checkbox"/>
---------------------------	-------	--------	-------	-----	--------------------------

Install the **Pressure Test Chamber** in the centrifuge.

Position 1:	_____	Pixels	(62 < value < 74)	OK:	<input type="checkbox"/>
Position 2:	_____	Pixels	Δ of position 1-2 _____ (390 < value < 410)	Pos 1-2 OK:	<input type="checkbox"/>
Position 3:	_____	Pixels	Δ of position 2-3 _____ (390 < value < 410)	Pos 2-3 OK:	<input type="checkbox"/>

6. Optical line sensor tests

Switch ON the SEPAX 2 and **wait 10 min.** before measuring


Attention: Calibration has to be done with the **lid closed** on the optical sensor.

Without filter

Red LED – Transmission	_____	V	(5.5 < value < 6,5V)	Red WF OK:	<input type="checkbox"/>
Blue LED – Transmission	_____	V	(6.1 < value < 7,1V)	Blue WF OK:	<input type="checkbox"/>

Black filter

Red LED – Transmission	_____	V	(0.0 < value < 0,2V)	Red BF OK:	<input type="checkbox"/>
Blue LED – Transmission	_____	V	(0.0 < value < 0,2V)	Blue BF OK:	<input type="checkbox"/>

	Work Procedure Maintenance	IT-1679-10
Created by: SAV		Page 3 / 4 21.01.2014
Approved: AQ		

7. Security tests

Check the splash detector (Attention: check the reactivity of the spill) approximately 2 sec	OK:	<input type="checkbox"/>
Check the three chamber cover sensors	OK:	<input type="checkbox"/>
Check the three temperature sensors	<div> <div>uC (20 ≤ T ≤ 40)</div> <div>Internal (20 ≤ T ≤ 40)</div> <div>Cover (15 ≤ T ≤ 35)</div> </div>	OK: <input type="checkbox"/>
For the cover temperature sensor, pass the finger on the reading window. The value must change	OK:	<input type="checkbox"/>
Main fan turns	OK:	<input type="checkbox"/>
There must be three beeps at start-up + red light	OK:	<input type="checkbox"/>

8. Centrifuge tower tests

Centrifuge tower alignment . The centrifuge should be perfectly centred	OK:	<input type="checkbox"/>
Centrifuge feet stability	OK:	<input type="checkbox"/>
Centrifuge upper deck (spill proof ring) is not touching the housing.	OK:	<input type="checkbox"/>
The gap must be minimum 0.7mm	OK:	<input type="checkbox"/>
Check if the centrifuge covers are neither broken nor bent, and do not create any vibrations	OK:	<input type="checkbox"/>
Run the centrifuge at 8000rpm with the pressure test chamber inside the cabinet, and check for vibrations.	OK:	<input type="checkbox"/>

9. Stepper motor tests

Stopcock configuration:	Standard <input type="checkbox"/>	Reduced (CS-570) <input type="checkbox"/>	
Run the protocol " test stopck pos " with stop cock ramp and choose time 60s			OK: <input type="checkbox"/>
Go to MAN mode:			OK: <input type="checkbox"/>
With the torque wrench , control the two stopcocks in both directions (0,3 < value < 0,5 Nm)			OK: <input type="checkbox"/>

Electrical test according to IEC 62353

Connect the Sepax2 to the Secutest using the power cable, which **belongs to the machine**. Hook up the Secutest to the sector (230VAC/110VAC).

10. Electrical earth test

In RCP mode

With the Secutest check the earth connection at 200mA. **Measuring points have to be < 0.3 Ω**

Test the following points:

The screw on the left hand side of the power connection	OK:	<input type="checkbox"/>
The rear fan grill	OK:	<input type="checkbox"/>
The top bag holder support (with bag holder inserted)	OK:	<input type="checkbox"/>
One central screw on the upper part of the head display	OK:	<input type="checkbox"/>
One screw underneath the Sepax (chassis)	OK:	<input type="checkbox"/>

11. Insulation test

In R-ISO mode

Set the Secutest to 500V. **Measuring points have to be above 70MΩ**


Test the following points:

The covers central knob	OK:	<input type="checkbox"/>
The luer connection	OK:	<input type="checkbox"/>
One of the four USB ports	OK:	<input type="checkbox"/>
One of the two Ethernet port	OK:	<input type="checkbox"/>
The handles fixing bolt on the left side of the device: Rear	OK:	<input type="checkbox"/>
The handles fixing bolt on the left side of the device: Front	OK:	<input type="checkbox"/>

12. Earth leakage test

In ICD mode

Current Normal Condition / Differential	_____	Max value < 0.5mA	OK: <input type="checkbox"/>
Single Fault Current / Alternatif	_____	Max value < 1 mA	OK: <input type="checkbox"/>

	Work Procedure Maintenance	IT-1679-10
		Page 4 / 4
Created by: SAV	Approved: AQ	21.01.2014

13. Current and voltage tests

In function mode

Current Consumption _____ Max value <0.65A @230Vac
 Max value <1A @110Vac OK: ☐

Other tests

14. Protocol tests

Put the value of position 1 in Section 5 at GMAP/Service menu/ Position piston offset. OK: ☐
 Run a complete procedure (e.g. UCB-HES) using a **water kit test**. OK: ☐

15. Visual control

Check the **exterior appearance** of the SEPAX 2 OK: ☐

16. AS-610 Traceability

Installed YES: ☐ NO: ☐
 Test the functionality of the barcode reader by scanning a barcode sample below OK: ☐ N/A: ☐
 Check the printer power supply and power cable exterior appearance OK: ☐ N/A: ☐
 Test the SEPAX 2 printing OK: ☐ N/A: ☐



17. Comments

18. Parts exchanged

Article Number	Designation	Old SN	New SN
In case of battery exchange		Reference	Manufacturing date

Tests conclusion

Tests with Automated Device Inspection Protocol passed PASS: ☐
 NOT PASS: ☐
 N/A: ☐

Tests passed for the SEPAX 2 PASS: ☐
 NOT PASS: ☐

Approval

Date: _____ Signed: _____